

Sandals & Heel Huggies



Size	Standard	Narrow
------	----------	--------

0	SL005	
1	SL007	SL017
2	SL009	SL019
3	SL011	SL021
4	SL013	SL023

****Narrow NOT available in size 0****

Cocoon Sandals Velcro



Size	Standard	Narrow
------	----------	--------

0	SL006	
1	SL008	SL018
2	SL010	SL020
3	SL012	SL022
4	SL014	SL024

****Narrow NOT available in size 0****

Cocoon Sandals Ski Buckle



Size	Standard	Narrow
------	----------	--------

1	N/A	SL045
2	N/A	SL047
3	N/A	SL049
4	N/A	SL051

****Can only be fitted with Velcro****

Cocoa Sandals Velcro



Size	Standard	Narrow
------	----------	--------

1	N/A	SL001
2	N/A	SL002
3	N/A	SL003
4	N/A	SL004

****Can only be fitted with Velcro****

Ferndale Sandals



Code SL053

This is a universal wedge used to fit all size sandals. They are stepped in 6mm, 12mm, 18mm & 24mm increments

Can be used with raisers

Depending on step used, longer bolts are supplied to attach them

Sandal Steps Wedge



Size	Standard	Narrow
------	----------	--------

1	SL025	SL030
2	SL029	SL034

****These are removable so you can get to the screws underneath****

Soft Insoles (Pair)

Made from Ninian vinyl & thick plastic to keep shape and strengthen

****All require 2 x M6 x 30 furniture bolts**

Sandals & Heel Huggies Cont.



Size	Standard	Narrow
0	SL035	SL040
1	SL036	SL041
2	SL037	SL042
3	SL038	SL043
4	SL039	SL044

****Can ONLY be fitted to cocoon sandals and CANNOT be retro fitted as they are stitched into the sandals at the time of ordering****

Sheepskin Cocoon Liner (Pair)



Size	12mm	18mm	24mm
0	SL055	SL060	SL065
1	SL056	SL061	SL066
2	SL057	SL062	SL067
3	SL058	SL063	SL068
4	SL059	SL064	SL069

****Will fit ALL sandals and a footplate is NOT required****

**** Please also note, there are also hundreds of RA codes for example RA045. SL codes are sandal raisers & RA codes are raisers used on footboards and ARE different products** ****

Sandal Raisers



Size	Code
1	HH001
2	HH002
3	HH003
4	HH004
5	HH005






****Velcro option as standard. If buckle is required, this can be done but at an additional cost to the customer.**

****If fitting heel huggie to ANY chair that's not fitted with a retractable base, a foot raiser IS required (RA113 & RA118)**


Heel Huggies

****All require 2 x M6 x 30 furniture bolts**


COCOON SANDALS (STD)




COCOON NARROW SANDALS




size 3 narrow



size 4 narrow




size 1 narrow




size 2 narrow

COCOA SANDALS


COCOA 3




COCOA 4








COCOA 1



COCOA 2



FERNDALE SANDALS



Heel Huggies Fitting Measurements

Product Codes	Depth	Max Circum	Min Circum
SIZE 1 HH001	50mm	160mm	130mm
SIZE 2 HH002	50mm	200mm	170mm
SIZE 3 HH003	75mm	240mm	210mm
SIZE 4 HH004	75mm	280mm	250mm
SIZE 5 HH005	90mm	310mm	280mm

Sandal Sizes

Sandal Model	Ferndale				Cocoa				Cocoon Narrow				Cocoon Standard				
Sandal size	1	2	3	4	1	2	3	4	1	2	3	4	0	1	2	3	4
Code	SL001	SL002	SL003	SL004	SL045	SL047	SL049	SL051	SL017	SL019	SL021	SL023	SL005	SL007	SL009	SL011	SL013
Shoe size (UK)																	
4	Available													Available			
5	Available				Available				Available					Available			
6	Available				Available				Available					Available			
7	Available	Available			Available					Available				Available			
8		Available			Available					Available				Available			
9		Available			Available						Available			Available	Available		
10		Available			Available	Available					Available			Available	Available		
11		Available	Available		Available	Available					Available			Available	Available		
12			Available		Available	Available					Available			Available	Available		
13			Available		Available	Available	Available				Available			Available	Available	Available	
1			Available		Available		Available				Available			Available		Available	
2			Available	Available	Available		Available				Available	Available		Available		Available	Available
3				Available	Available		Available				Available	Available		Available		Available	Available