

juni



our vibrant chair providing unique contouring, for mild to moderate postural support from 18mths old

**Numerous benefits**

Whether in the classroom, nursery, or at home, better posture whilst seated creates a more positive experience for your child; fine motor skills develop, confidence increases and the ability to engage and learn becomes easier. The Juni comes in three vibrant colours and the contemporary appeal helps to promote social inclusion.

**Posture at its foundation**

Juni's seat and backrest both provide unique contouring. This promotes better posture, gives higher levels of comfort and helps your child feel more secure and stable in the Juni. It is our vibrant chair for children who require minimal levels of postural support when seated.

**Highly adjustable**

We understand that your child's needs change during each stage of growth and development. Seat height and depth as well as armrest height can all be easily adjusted. The backrest can be reclined by 5° increments (from 95° to 105°), whilst it also features an unique built-in lumbar support.

high levels of comfort



**Standard key features**



Angle adjustable back rest



Adjustable lumbar support



Easy leg adjustments



Lap strap



**Juni technical data**

all dimensions in mm unless otherwise stated

| Size                 | 1          | 2          | 3          | 4          |
|----------------------|------------|------------|------------|------------|
| Height               | 160-220    | 200-280    | 240-340    | 300-400    |
| Width                | 240        | 280        | 325        | 365        |
| Depth                | 180-240    | 220-300    | 260-360    | 320-420    |
| Back height          | 260        | 290        | 325        | 400        |
| Back angle (°)       | 95/100/105 | 95/100/105 | 95/100/105 | 95/100/105 |
| Max user weight (kg) | 30         | 40         | 50         | 60         |
| Product weight (kg)  | 4          | 6          | 8          | 10         |

| Colour options | Codes |       |       |       |
|----------------|-------|-------|-------|-------|
| Red            | 16011 | 16012 | 16013 | 16014 |
| Lime           | 16021 | 16022 | 16023 | 16024 |
| Blue           | 16031 | 16032 | 16033 | 16034 |
| Black          | 16041 | 16042 | 16043 | 16044 |

The Juni chair comes with contoured seat and backrest, built-in lumbar support, pelvic lap strap, height and depth adjustable seat, height adjustable armrests as well as angle adjustable backrest as standard.

Instructional video available online

# juni

## optional accessories



### Skis

CODE SK101 - SK104

Helps to stabilise and prevent chair movements



### Ski & footboard<sup>1</sup>

CODE SK111 - SK114

Provides foot support and offers greater stability



### Mobile skis

CODE DL101 - DL104


Gives clearance to facilitate standing transfer



### Sliding footrest<sup>2</sup>

CODE DL111 - DL114

Provides foot support  
Slides backwards to give clearance to facilitate standing transfer



### Quick release tray

Patent applied for  
CODE TR101 - TR104

The tray is available in white



### Height adjustable push handle

CODE PH101 - PH104

Adjusts to achieve optimum push position




### Side support cushions<sup>3</sup>

Provides hip and thigh support  
Available in 25mm and 50mm



### Footrest

CODE FR051 - FR054  
Available in 4 sizes



### 4-point harness<sup>3</sup>

Provides trunk support  
Available in 3 colours & 4 sizes



### Pommel

Code PL111 - PL132  
Encourages hip abduction

<sup>1</sup> Raisers are also available in either 12mm [CODE: RA049 - RA052] or 24mm [CODE: RA045 - RA048]

<sup>2</sup> Sliding footrest can only be used in conjunction with the mobile skis

<sup>3</sup> Codes vary depending on the colour and size; please contact us directly or see our price list for all codes

# juni

## all optional accessories



|   | Size 1<br>Code | Size 2<br>Code | Size 3<br>Code | Size 4<br>Code |
|---|----------------|----------------|----------------|----------------|
| <b>Lower leg &amp; seat accessories</b>     |                |                |                |                |
| skis  | SK101          | SK102          | SK103          | SK104          |
| mobile skis                                 | DL101          | DL102          | DL103          | DL104          |
| skis & footboard                            | SK111          | SK112          | SK113          | SK114          |
| sliding footrest (for use with mobile skis) | DL111          | DL112          | DL113          | DL114          |
| height adjustable push handle               | PH101          | PH102          | PH103          | PH104          |
| juni footrest                               | FR051          | FR052          | FR053          | FR054          |
| juni extended outriggers                    | LG153          | LG154          | LG155          |                |
| pommel - black                              | PL141          | PL141          | PL141          | PL142          |
| <b>Tray accessories</b>                     |                |                |                |                |
| quick release tray                          | TR101          | TR102          | TR103          | TR104          |
| <b>Cushions</b>                             |                |                |                |                |
| 25mm side support cushions - red            | CS511          | CS511          | CS512          | CS512          |
| 25mm side support cushions - lime           | CS521          | CS521          | CS522          | CS522          |
| 25mm side support cushions - blue           | CS531          | CS531          | CS532          | CS532          |
| 25mm side support cushions - black          | CS541          | CS541          | CS542          | CS542          |
| 50mm side support cushions - red            | CS611          | CS611          | CS612          | CS612          |
| 50mm side support cushions - lime           | CS621          | CS621          | CS622          | CS622          |
| 50mm side support cushions - blue           | CS631          | CS631          | CS632          | CS632          |
| 50mm side support cushions - black          | CS641          | CS641          | CS642          | CS642          |
| 25mm + 25mm double cushions - red           | JC001          | JC001          | JC004          | JC004          |
| 25mm + 25mm double cushions - lime          | JC002          | JC002          | JC005          | JC005          |
| 25mm + 25mm double cushions - blue          | JC003          | JC003          | JC006          | JC006          |
| 25mm + 25mm double cushions - black         | JC007          | JC007          | JC008          | JC008          |
| 50mm + 25mm double cushions - red           | JC051          | JC051          | JC054          | JC054          |
| 50mm + 25mm double cushions - lime          | JC052          | JC052          | JC055          | JC055          |
| 50mm + 25mm double cushions - blue          | JC053          | JC053          | JC056          | JC056          |
| 50mm + 25mm double cushions - black         | JC057          | JC057          | JC058          | JC058          |
| <b>Harnesses</b>                            |                |                |                |                |
| 4 point harness - red                       | HN211          | HN212          | HN213          | HN214          |
| 4 point harness - lime                      | HN221          | HN222          | HN223          | HN224          |
| 4 point harness - blue                      | HN231          | HN232          | HN233          | HN234          |
| 4 point harness - black                     | HN251          | HN252          | HN253          | HN254          |

### Dolly Base

CODE PH101 - PH104

Facilitates manoeuvrability of your Juni chair

Available in 4 sizes



Smirthwaite's mobile base for your Juni, opens up a whole new world.

- Manoeuvre from class to class or around the home.
- Slimline and easy to use.
- Hard wearing wheels.
- Compliments the style of your Juni.

Colour Ranges  
See page 55 for colour options

# juni extras

We work in close collaboration with therapists, parents and children to ensure that our products meet and surpass your expectations. Through your valuable feedback the Juni also has the following optional extras: highly adjustable Juni Table, Juni Hi & Lo Boosters and Juni Hi-lo Base to overcome height issues at home, nursery & school. See all Juni accessories on pages 38 and 39.

## Juni Hi-lo Base

We understand the challenges of the classroom and home environment. Height adjustment is now easy with our Juni Hi-lo Base using the gas assisted lift. The slimline footprint of the base makes this the ideal solution to fit under most classroom desks. So now, you can use your Juni chair together with the Juni Hi-lo Base by simply attaching it to the plate on top of the base; or use the chair independently where height adjustments are not needed.

Only suitable for Juni chairs size 1 & 2.



## Juni Hi-lo Base technical data

all dimensions in mm unless otherwise stated

| Size of Juni                 | 1          | 2          |
|------------------------------|------------|------------|
| Seat height from floor       | 160 - 570  | 200 - 570  |
| Lower leg length to footrest | 160 - 220  | 200 - 280  |
| Push handle height           | 500 - 1000 | 500 - 1000 |
| Base length and width        | 650 x 520  | 650 x 520  |
| Max user weight (kg)         | 30kg       | 40kg       |
| Product weight (kg)          | 10.5kg     | 10.5kg     |

| Colour           | 1     | 2     |
|------------------|-------|-------|
| Red footrest     | HL111 | HL112 |
| Lime footrest    | HL121 | HL122 |
| Blue footrest    | HL131 | HL132 |
| Natural footrest | HL141 | HL142 |
| Black footrest   | HL151 | HL152 |

Colour Ranges  
See page 55 for colour options

Juni extras to help overcome height issues at home, nursery & school

# juni extras

## Juni Table

The Juni table offers two key benefits; the legs are adjustable so your child can work or play at the correct height and the table folds, which allows for easy storage. Available in a variety of colours this is the perfect companion to your favourite Juni chair.



## Juni table technical data

all dimensions in mm unless otherwise stated

| Size                | 1       | 2       |
|---------------------|---------|---------|
| Height              | 310-450 | 400-600 |
| Width               | 700     | 780     |
| Depth               | 530     | 590     |
| Max weight (kg)     | 60      | 60      |
| Product weight (kg) | 8       | 10      |

| Colour      | Code |      |
|-------------|------|------|
| Red top     | 2411 | 2412 |
| Lime top    | 2421 | 2422 |
| Blue top    | 2431 | 2432 |
| Natural top | 2441 | 2442 |



## Juni Boosters

Our Hi Booster is perfect for the breakfast bar or for feeding your child without crouching or kneeling. Our Lo Booster allows for sitting your child at the dinner table or giving your child extra height.

Hi & Lo Boosters come complete with a black 5-point harness and are also available with castors\*. The Juni Boosters accommodate Juni chair size 1 and 2 only. Please visit our websites for more details.

## Juni Boosters technical data

all dimensions in mm unless otherwise stated

| Booster Type             | Lo Booster      | Hi Booster      |
|--------------------------|-----------------|-----------------|
| Size of chair on booster | 1 & 2           | 1 & 2           |
| Seat height              | 670-770         | 870-970         |
| Footprint open           | 520 x 520       | 520 x 640       |
| Folded size              | 520 x 340 x 700 | 520 x 420 x 930 |
| Max user weight (kg)     | 40              | 40              |
| Product weight (kg)      | 7               | 8               |

| Colour  | Size   |        | Codes  |        |
|---------|--------|--------|--------|--------|
|         | Size 1 | Size 2 | Size 1 | Size 2 |
| Red     | JB111  | JB121  | JB112  | JB122  |
| Lime    | JB211  | JB221  | JB212  | JB222  |
| Blue    | JB311  | JB321  | JB312  | JB322  |
| Natural | JB411  | JB421  | JB412  | JB422  |
| Black   | JB511  | JB521  | JB512  | JB522  |

