

# whirl packages

a popular prone stander with a unique module design, promoting excellent symmetry

## Supportive

The improved Whirl enables children from 12 months old to young adults to gradually develop and maintain their ability to weight bear by promoting excellent symmetry. It is suitable for children who require lower to higher levels of postural support and control. The Whirl is ideal for both school and home environments.


## Packages

All options come complete with:

- anterior, pelvic and knee supports
- cocoon sandals,
- tray grab rail full width
- activity tray
- bowl with cover



## Deluxe Colour Range



## Whirl technical data

all dimensions in mm unless stated

| Size   | 1         | 2          |
|--|-----------|------------|
| Code   | 8353      | 8354       |
| Thoracic support (top) to foot plate height    | *630-1020 | *760-1270  |
| Pelvic support (top) to foot plate height      | *410-830  | *560-1030  |
| Knee support (top) to foot plate height        | *130-420  | *165-520   |
| Foot plate to floor height (excluding sandals) | 70        | 70         |
| Max width (between thoracic supports)          | 300       | 460        |
| Tray size (depth and width)                    | 610 x 690 | 700 x 800  |
| Angle (°)                                      | 0 to 30   | 0 to 30    |
| Base frame                                     | 890 x 660 | 1020 x 710 |
| Max user weight (kg)                           | 70        | 90         |
| Product weight (kg)                            | 38        | 40         |

\*Minimum heights taken with footplate raised 50mm from floor

Instructional video available online

+44 (0)1626 835552  
info@smirthwaite.co.uk  
www.smirthwaite.co.uk

Supporting children &  
their families for life



**Smirthwaite Ltd**  
16 Wentworth Road, Heathfield, Newton Abbot, Devon TQ12 6TL