

★ mobi heavy  
duty changer

**1** designed to meet the needs of plus sized users, includes an powerful motor and extra wide mattress



**MADE IN BRITAIN**



[smirthwaite.co.uk](http://smirthwaite.co.uk)

★ mobi heavy  
duty changer



## key features

- ★ Electronically operated
- ★ Travel range of 370mm
- ★ Weight capacity up to 500kg
- ★ Mobile & easy to manoeuvre
- ★ Stowaway side rails
- ★ Dry room environments only



## key benefits

- ★ **An adult changing table** with larger weight requirements. Its powerful 12K Newton motor, coupled with the extra wide mattress platform, makes it the ideal solution.
- ★ **Does not require any wall mounting or floor fixing.** They are easy to manoeuvre which means the user can be transported from room to room with no trouble (door opening has to be 1100mm or more).
- ★ **Electronically operated**, which makes adjusting the height of the table simple, this is beneficial for both the carer and user.
- ★ **Levitex technology** is available as an added extra.

☎ 44 (0)1626 835552

✉ [info@smirthwaite.co.uk](mailto:info@smirthwaite.co.uk)

