


# harness measurements

## heathfield, brookfield and combi harnesses



**Teenage harness**

Size	1	3
A (cm)	31	36
B (cm)	28	34
C (cm)	18	23
D (cm)	15	16



**Zipped jacket**

Size	1	2	3
A (cm)	7	9	11
B (cm)	24	29	33
C (cm)	25.5	32	35.5
D (cm)	7	9	11
E (cm)	5	5	5.5




**Butterfly harness**

Size	One Size
A (cm)	15
B (cm)	13



**Bib harness (padded)**


Size	1	2	3
A (cm)	11	14	17.5
B (cm)	23	26.5	32
C (cm)	22	28	31
D (cm)	9	13	15.5
E (cm)	6	6.5	7.5



**Waistcoat harness**


Size	1	2	3
A (cm)	5	10	13
B (cm)	15	20	25.5
C (cm)	31	39	40
D (cm)	10	15	22
E (cm)	5	5	6.5
F (cm)	63	78	98

F = circumference front to back



**Groin harness**

Size	1	2	3
A (cm)	24	27	32
B (cm)	17.5	25	28
C (cm)	37.5	47	60
D (cm)	6.5	8	9
Weight	150g	200g	300g




**4-point harness**

Size	1	2	3
A (cm)	23	28	29
B (cm)	28	28	30

# harness measurements




## samba and hardrock harnesses



**Soft cross harness**

Size	1	2	3	4
A (cm)	46	51	57	63
B (cm)	5	5	5	5



**Shoulder harness for the mid-line harness**


Size	1	2	3	4
A (cm)	22	22	22	28
B (cm)	16	16	16	22
C (cm)	16	16	16	16
D (cm)	5	5	5	7
E (cm)	27	27	27	27
F (cm)	5	5	5	6



**Mid-line support harness with sheepskin pad**

Size	1	2	3	4
A (cm)	16	16	16	20
B (cm)	12	12	12	12
C (cm)	7 - 17	7 - 17	7 - 17	7 - 17

## indigo harnesses



**Teenage harness  
(with or without zip)**

Size	1
A (cm)	36
B (cm)	34
C (cm)	23
D (cm)	16

Smirthwaite Ltd  
16 Wentworth Road  
Heathfield Industrial Estate  
Newton Abbot  
Devon  
TQ12 6TL

[www.smirthwaite.co.uk](http://www.smirthwaite.co.uk)  
[info@smirthwaite.co.uk](mailto:info@smirthwaite.co.uk)

T: 01626 835552

Smirthwaite harnesses offer comfort whilst maintaining secure positioning in all our seating systems. A range of colours and fabrics are available and 'made to measure' solutions and modifications can be made upon request.

Designed for use on Smirthwaite products only